

BD Microtainer™ Plasma Separator Tube

Amber reservoir with Lithium Heparin
and Microgard™ Closure



BD

Indispensable to
human health

Designed For Light Sensitive Tests

BD Microgard™ Closure

Recessed Sealing Plug

Note: The BD Microgard™ Closure has been cut away to demonstrate its unique features.

Integral Collector

Wider Tube Diameter


Textured Inner Surface

Printed Lot #'s, Fill Lines
and Expiration Date

Precisely Measured
Lithium Heparin

Inert Polymer Gel

Amber colored tube for
light sensitive analyte testing
(i.e. bilirubin tests)



Designed for convenience...

- Integral collector eliminates the need for assembly
- Lot number and expiration dates printed on tubes, for increased quality control
- 8-mm inner tube diameter allows for easy probe access
- BD Microgard™ Closure nests securely on bottom of tube or extender

performance...

- Wider tube reservoir with exclusive textured inner surface improves blood flow and mixability
- Precisely measured and dispensed lithium heparin

... and safety

- Unique "twist-assist" safety BD Microgard™ closure significantly reduces splatter upon removal

BD Microtainer™ Plasma Separator Tube

Amber reservoir with Lithium Heparin
and Microgard™ Closure

ORDERING INFORMATION

	Product Description	Reorder Number	Quantity per case	Quantity per bag
	Amber BD Microtainer™ Plasma Separator Tube with Lithium Heparin and Microgard™ Closure	365987	200	50
	Microtainer™ Tube with Lithium Heparin and Microgard™ Closure	365965	200	50
	Microtainer™ Tube with EDTA and Microgard™ Closure	365974	200	50
	Microtainer™ Tube Extender	365976	200	50

Caution:

Handle all biologic samples and blood collection products in accordance with the policies and procedures of your facility.

Obtain appropriate medical attention in the event of any exposure to biologic samples (for example, through a puncture injury) since the samples may transmit HBV (hepatitis), HIV (AIDS), or other infectious diseases.

Discard all products contaminated with blood in biohazard containers approved for their disposal.

Packaging designed to store more tubes in less space
For more information, please call customer service
1-888-237-2762



BD Vacutainer Systems
Preanalytical Solutions
1 Becton Drive
Franklin Lakes, NJ 07471
www.bd.com